


Welcome to Barcoding the Broads

A Wellcome-funded programme of public engagement events and activities to explore biodiversity on the Norfolk Broads, led by the Earlham Institute as part of our work on the Darwin Tree of Life project.

Content:

- 1 DNA Barcoding Protocol - Invertebrates - All Equipment, Consumables and Reagents**
- 2 Stage 1] Sample Collection and DNA Extraction - Equipment list**
- 3 Stage 1] Sample Collection and DNA Extraction**
- 4 Stage 2] Polymerase Chain Reaction (PCR) - Equipment list**
- 5 Stage 2] Polymerase Chain Reaction (PCR)**
- 6 Stage 3] Gel Electrophoresis and Sanger Sequencing - Equipment list**
- 7 Stage 3] Gel Electrophoresis and Sanger Sequencing**
- 8 Stage 4] - DNA Subway**
- 9 Metadata Collection Table**



DNA Barcoding Protocol - Invertebrates



All Equipment, Consumables, and Reagents

| Equipment: | Check |
|---|-------|
| -20 °C freezer | |
| Balance | |
| BlueGel electrophoresis system | |
| Centrifuge | |
| Computer/Laptop with internet connection | |
| Conical flask | |
| Eppendorf and PCR tube rack | |
| Ice bucket | |
| Lab book | |
| Media bottle | |
| Metal forceps | |
| Metal spatula | |
| Microwave | |
| MiniPCR thermal cycler | |
| Nitrile gloves | |
| 10 µL, 20 µL, 200 µL and 1000 µL pipettes | |
| Plastic pestles | |
| Scissors | |
| Waste bin | |
| Water bath (or mug and kettle) | |

| Consumables: | Check |
|---|-------|
| 0.2 mL PCR tubes | |
| 1.5 mL eppendorf tubes | |
| 10 mL plastic syringe | |
| 10 µL, 20 µL, 200 µL and 1000 µL pipette tips | |
| Aluminium foil | |
| Ice | |
| Marker pen | |
| Parafilm M | |
| Plastic weigh boats | |
| Sanger sequencing labels (Genewiz/Azenta) | |
| Sterile/Distilled water | |

| Reagents: | Check |
|--|-------|
| 1 x TBE buffer [NB: if using 10 x TBE buffer, dilute ten-fold in sterile/distilled water] | |
| 10% chelex solution [made by mixing 1g chelex into 9 mL sterile/distilled water] | |
| 100bp DNA ladder | |
| Agarose | |
| DNA loading dye | |
| Ethanol (for cleaning) | |
| CO1 primers (forward and reverse) | |
| Ready-to-Go PCR beads | |
| SYBR Safe DNA gel stain | |

DNA Barcoding Protocol - Invertebrates - STAGE 1



Sample Collection and DNA Extraction



| Equipment: | Check |
|--------------------------------|-------|
| Balance | |
| Centrifuge | |
| Eppendorf and PCR tube rack | |
| Lab book | |
| Media bottle | |
| Metal forceps | |
| Metal spatula | |
| Nitrile gloves | |
| 200 μ L pipette | |
| Plastic pestle | |
| Scissors | |
| Waste bin | |
| Water bath (or mug and kettle) | |

| Reagents: | Check |
|---|-------|
| 10% chelex solution [made by mixing 1g chelex into 9mL sterile/distilled water] | |
| Ethanol (for cleaning) | |

Notes



| Consumables: | Check |
|--------------------------|-------|
| 1.5 mL eppendorf tubes | |
| 10 mL plastic syringe | |
| 200 μ L pipette tips | |
| Aluminium foil | |
| Marker pen | |
| Parafilm M | |
| Plastic weigh boats | |
| Sterile/Distilled water | |

DNA Barcoding Protocol - Invertebrates - STAGE 1



Sample Collection and DNA Extraction



The first task is to collect and document an invertebrate sample, ensuring you respect the study environment and have permission to collect in the chosen location. Only a small amount of material is needed for DNA barcoding.

Metadata



As much information as possible about the sample should be recorded alongside the collection.

This includes:

- Photos of the invertebrate
- Latitude, longitude, temperature and elevation of the study location
- Time and date of collection
- Description of the invertebrate and the habitat it was found
- Sample number and name of the collector

All sample information should be logged in a lab book or on a computer.

An example metadata collection table can be found on **page 9**



DNA Extraction



Wear gloves throughout to avoid contamination



Use a fresh tip every time you pipette to avoid contamination



DNA is extracted from the invertebrate and separated from contaminants using physical and chemical methods.

- 1** Add 100 μL of 10 % chelex solution to a 1.5 mL eppendorf tube and use a marker pen to label the top and side of the tube with the sample name/number.
- 2** Carefully remove a small section of your invertebrate sample, wearing gloves to avoid contamination - e.g. smaller specimens like ants and spiders can be halved or used whole, or a few millimeters of leg/antenna can be taken from larger specimens like beetles and butterflies
- 3** Transfer the section of invertebrate sample to the labelled eppendorf tube and crush thoroughly into the chelex solution using a plastic pestle until it is completely broken apart
- 4** Close the eppendorf tube and seal it tightly with parafilm. Heat the contents for 10-15 mins in a water bath (85-90 $^{\circ}\text{C}$) or in a mug containing boiling water. Avoid contamination by ensuring the lid of the tube isn't submerged.
- 5** Use a centrifuge (e.g. 10000 rpm for 2 min) to pellet the solid material. Alternatively, allow the contents to settle for 15 mins until two distinct layers form. The DNA is found in the top layer of liquid that forms above the solid material. Handle the tube carefully to avoid disturbing the solid layer at the bottom

DNA Barcoding Protocol - Invertebrates - STAGE 2



Polymerase Chain Reaction (PCR)



| Equipment: | Check |
|--------------------------------------|-------|
| -20 °C freezer | |
| Computer/Laptop for miniPCR software | |
| Eppendorf and PCR tube rack | |
| Ice bucket | |
| Lab book | |
| MiniPCR thermal cycler | |
| Nitrile gloves | |
| 20 µL and 200 µL pipettes | |
| Waste bin | |

| Consumables: | Check |
|-------------------------------|-------|
| 0.2 mL PCR tubes | |
| 20 µL and 200 µL pipette tips | |
| Ice | |
| Marker pen | |
| Plastic weigh boats | |
| Sterile/Distilled water | |

| Reagents: | Check |
|-----------------------------------|-------|
| Ready-to-Go PCR beads | |
| C01 primers (forward and reverse) | |
| Ethanol (for cleaning) | |

Notes



DNA Barcoding Protocol - Invertebrates - STAGE 2



Polymerase Chain Reaction (PCR)



Wear gloves throughout to avoid contamination



Use a fresh tip every time you pipette to avoid contamination



PCR is used to amplify (make billions of copies of) a specific region of DNA. For invertebrates, the DNA barcode region is called **CO1**.

- 1** Transfer a ready-to-go PCR bead (containing Taq polymerase, nucleotides, stabilisers and reaction buffer) to a 0.2 mL PCR tube and use a marker pen to label the top and side of the tube with the sample name/number
- 2** Add 20 μ L sterile/distilled water to the PCR tube and allow the bead to completely dissolve
- 3** Add 1.5 μ L of the forward primer and 1.5 μ L of the reverse primer to the PCR tube. Primer sequences are:
CO1-M13-forward: 5'-TGTAACGACGGCCAGTGGTCAACAAATCATAAAGATATTGG-3'
CO1-M13-reverse: 5'-CAGGAAACAGCTATGACTAACTTCAGGGTGACCAAAAAATCA-3'
- 4** Add 2 μ L of extracted DNA (from the top layer of liquid that forms above the solid material at the end of STAGE 1) to the PCR tube. Mix thoroughly with a pipette tip before firmly closing the lid. The final volume of the solution in the PCR tube is 25 μ L. A negative control (without DNA) is prepared by repeating the steps above using 2 μ L sterile/distilled water instead of 2 μ L extracted DNA
- 5** Transfer the PCR tubes to the miniPCR machine. Screw the lid of the machine closed and then start the thermal cycling protocol using the miniPCR software

Thermal Cycling



Amplification of the DNA barcode region happens in three steps:

- 1** *Denaturation*, where the double-stranded DNA is heated to separate the strands
- 2** *Annealing*, where the primers mark/flank the DNA barcode region by binding to the separated strands of DNA
- 3** *Extension*, where the Taq polymerase synthesises new DNA strands.

These steps are repeated (cycled) to exponentially produce exact copies of the DNA barcode region.

Thermal cycling programme for CO1 (1 hr 57 min duration, with heated lid) is:

Initial denaturation = 94 °C for 60 sec

| *Denaturation* = 95 °C for 30 sec

| *Annealing* = 50 °C for 30 sec

| *Extension* = 72 °C for 45 sec

| Number of cycles = 35

Final extension = 72 °C for 180 sec

DNA Barcoding Protocol - Invertebrates - STAGE 3



Gel Electrophoresis and Sanger Sequencing



| Equipment: | Check |
|---|-------|
| -20 °C freezer | |
| 10 µL, 20 µL, 200 µL and 1000 µL pipettes | |
| Balance | |
| BlueGel electrophoresis system | |
| Computer/Laptop with internet connection | |
| Conical flask | |
| Eppendorf and PCR tube rack | |
| Ice bucket | |
| Lab book | |
| Metal spatula | |
| Microwave | |
| Nitrile gloves | |
| Waste bin | |

| Consumables: | Check |
|---|-------|
| 0.2 mL PCR tubes | |
| 1.5 mL eppendorf tubes | |
| 10 mL plastic syringe | |
| 10 µL, 20 µL, 200 µL and 1000 µL pipette tips | |
| Ice | |
| Marker pen | |
| Parafilm M | |
| Plastic weigh boats | |
| Sanger sequencing labels (Genewiz/Azenta) | |
| Sterile/Distilled water | |

| Reagents: | Check |
|---|-------|
| 1 x TBE buffer [NB: if using 10 x TBE buffer , dilute ten-fold in sterile/distilled water] | |
| 100bp DNA ladder | |
| Agarose | |
| DNA loading dye | |
| Ethanol (for cleaning) | |
| SYBR Safe DNA gel stain | |

Notes



DNA Barcoding Protocol - Invertebrates - STAGE 3



Gel Electrophoresis



Use a fresh tip every time you pipette to avoid contamination



Wear gloves throughout to avoid contamination



Gel electrophoresis confirms whether the polymerase chain reaction (PCR) was successful. The technique mobilises DNA by charge and DNA fragments are separated by size, with smaller fragments moving further through the gel.

- 1** Mix 0.4 g agarose with 20 mL of 1 x Tris/Borate/EDTA (TBE) buffer in a conical flask. NB: if using **10 x TBE buffer**, remember to dilute ten-fold in sterile/distilled water
- 2** Microwave the mixture for 20-30 seconds or until the agarose has completely dissolved, ensuring the solution does not boil over. Once the solution has cooled slightly, add 2 µL of SYBR Safe DNA gel stain (NB: reagent is light sensitive so store in the dark)
- 3** While it's still warm, gently swirl the solution and carefully pour into the gel tray with a comb insert, giving 9 large wells (20 µL capacity) or 13 small wells (10 µL capacity). Allow the gel to set completely (approximately 15 minutes)
- 4** Once the gel has set, carefully remove the comb and transfer the gel tray into the main sample tray. Cover the gel with 35 mL of 1 x TBE buffer (NB: if using **10 x TBE buffer**, remember to dilute ten-fold in sterile/distilled water). DNA samples can now be prepared for analysis on the gel
- 5** Add 5 µL of amplified DNA sample (post-PCR) to 1 µL loading dye in a 0.2 mL PCR tube and mix with a pipette tip. Repeat with the negative control sample. The DNA ladder (used to determine the approximate size of DNA samples on the gel) is similarly prepared by combining 5 µL of 100bp DNA ladder with 1 µL loading dye in a 0.2 mL PCR tube
- 6** Carefully add the 6 µL of each DNA/dye sample to the wells on the gel (one sample per well) then seal the tank with the lid and press the power button. Progress of the gel is monitored by turning on the blue light (viewed through the hole in the dark box)

Sanger Sequencing



DNA (Sanger) sequencing is used to determine the order of A, C, T & G building blocks in the DNA barcode region. If gel electrophoresis confirms the PCR experiment was successful then 20 µL of the amplified DNA sample can be transferred to a 1.5 mL Eppendorf tube and logged for collection/analysis via the Genewiz/Azenta website.

Progress with sample analysis is tracked via the Genewiz/Azenta website and **.ab1** files for the DNA Subway can be downloaded once the process is complete. Speak to your local Genewiz/Azenta account manager for more information about submitting samples for DNA sequencing and reviewing the results.

DNA Barcoding Protocol - Invertebrates - STAGE 4



DNA Subway



Equipment



Computer/Laptop with internet connection

Software



DNA sequences are analysed using the "Blue Line" on the DNA Subway (<https://dnasubway.cyverse.org>). The platform is accessed with a CyVerse account or as a Guest (for the latter, projects can't be saved long-term). Select "Barcoding: CO1" as the project type and upload **.ab1** files from the Genewiz/Azenta website (one forward and one reverse per sample). A project title and short description can be added here too.

Step 1: Assemble Sequences

DNA sequences can be viewed and automatically trimmed before the forward and reverse sequences for each sample are paired up. A consensus sequence is then built from the forward/reverse pairs and can be manually trimmed at both ends, ideally keeping the final sequence length above 400bp.

Step 2: Add Sequences

The consensus sequence can be compared with online databases using a tool called BLAST (Basic Local Alignment Search Tool) to identify invertebrate species from the CO1 DNA barcode region. Suggested species names are shown alongside an analysis of how close the consensus sequence matches those in the online databases. A range of CO1 reference data sets (e.g. common insects) can be added to compare against the consensus sequence in the next step.

Step 3: Analyse Sequences

Multiple sequences can be compared using a tool called MUSCLE (Multiple Sequence Comparison by Log Expectation). Grey areas indicate a match with the overall consensus whereas coloured bars/letters show where individual sequences differ. The multiple sequence alignments can be trimmed and then phylogenetic trees can be made to visualise the relationship between invertebrate species. Two types of phylogenetic tree are available: neighbour joining (NJ) and maximum likelihood (ML). For both, an outgroup (i.e. the invertebrate species in the data set that is most dissimilar to the one being studied) is chosen to act as a reference point for determining where the root and branches of the tree are placed.

➡ More information about the DNA Subway can be found under the "DNA Sequencing and Analysis" tab on the Barcoding the Broads webpage. Please log your DNA barcoding results on the Darwin Tree of Life community map/database to support sampling for full genome analysis and request assistance with species identification.

Barcoding the Broads Metadata Collection Table



Name: _____ Date (dd/mm/yy): _____ Sample (plant/invertebrate/fungus): _____

Organisation/School: _____ Collection Location: _____

| Specimen Collection Data | | | | | | | | |
|--------------------------|--------------------------|------|----------|-----------|-----------|-------------------|---------------------|-------|
| Sample Number | Photo taken | Time | Latitude | Longitude | Temp (°C) | Predicted Species | Habitat Description | Notes |
| | <input type="checkbox"/> | | | | | | | |
| | <input type="checkbox"/> | | | | | | | |
| | <input type="checkbox"/> | | | | | | | |
| | <input type="checkbox"/> | | | | | | | |