




# Make up 4 boxes




Place a guidance sheet into each box




Take your strip of barcodes (x8) and cut these into x4 pairs, of which the 2D barcodes should match across the pair. Check participant ID matches across all four pairs.




Place one pair of barcode labels into each box




Place a safety bag inside each box



Place an absorbent sheet into a griplock bag and place inside the box (repeat for all boxes in your set of four)



Place a swab kit inside each box



Close all four kit boxes and place them in a single carrier bag. Double knot the handles and drop into your completed kit receptable.